



New and Used Medical Equipment For Sale and Wanted
We don't own it... We help you sell it.

[My DOTmed](#)
[Login](#) | [Register](#)

[All](#) [Auctions](#) [Equipment](#) [Parts](#) [Careers](#) [Services Directory](#) [News](#) [Forums](#)

The Medical Industry
Business Weekly

[?]

[Search Forums/News](#)

[DOTmed Home](#) [News Home](#) [Industry](#) [People & Companies](#) [Shows & Events](#) [Forums](#) [Career Corner & Jobs](#)

July 10, 2008

Other Headlines

[Senate Agrees to Restore Doctors' Pay](#)

Senator Kennedy's dramatic return results in key vote in the measure to return physician Medicare payment to pre-cut levels.

[Siemens Job Cuts Set at 16,750](#)

Company confirms it will cut about 4% of its workforce worldwide; healthcare division to lose 2,800 positions.

[DOTmed Certified Sharon Moore of Capital Medical Equipment Sees Expansion Down the Road](#)

After 30 years in a hospital setting, she has joined the ranks of imaging saleswomen.

[Florida Event a Major Happening](#)

Florida International Medical Exposition (FIME) 2008 will help exhibitors gain a competitive edge in the marketplace; Visit DOTmed at Booth #1511.

[Strategic Medical Sales -- Experience Makes the Difference](#)

Expertise and know-how attract new domestic and international customers for patient monitor sales.

Have News for Us?

[Submit your news](#) on the industry, people, or companies.

New MRI Developed by Research Team

by Joan Trombetti, Writer

Researchers at the Department of Energy's Lawrence Berkeley National Laboratory, and the University of California at Berkeley, have developed a magnetic resonance imaging method that is faster, more selective, and many thousands of times more sensitive than standard MRI.



Leif Schroeder (left) and team members

Using "temperature-controlled molecular depolarization gates," the new technique builds on a series of previous developments in MRI, and the closely related field of nuclear magnetic resonance (NMR).

Standard MRI suffers from low sensitivity and requires patients to remain motionless for long periods of time inside noisy, claustrophobic machines. The new method requires less time and is able to distinguish even among specific target molecules.

"The new method holds the promise of combining a set of proven NMR tools for the first time into a practical, supersensitive diagnostic system for imaging the distribution of specific molecules on such targets as tumors in human subjects or even on individual cancer cells," said lead author Leif Schroeder of Berkeley Lab's Materials Sciences Division.

MRI and NMR make use of the quantum-mechanical phenomenon known as nuclear spin; nuclei with odd numbers of protons or neutrons will orient themselves like tiny bar magnets, spin "up" or spin "down," in a strong magnetic field. If the spinning nuclei are knocked off-axis by a jolt of radio-frequency (RF) energy, they wobble or precess at a characteristic rate, a rate that is strongly conditioned by their immediate chemical neighbors. During a certain relaxation time (typical of each atomic species in a specific environment), the nuclei reorient themselves and emit a radio signal that reveals both their position and their chemical surroundings.

The spin-up state requires fractionally less energy, so there's typically a slight excess of spin-up nuclei, about one in a hundred thousand (.001 percent) and it's this tiny difference that yields a useful signal. In clinical settings MRI is usually done using hydrogen nuclei, protons, which are ubiquitous in the human body. But other nuclear species, notably the noble gas xenon, offer advantages over hydrogen that in the case of xenon include a virtual absence of background signal, since there is no xenon in biological systems.

Xenon is useful in MRI and NMR because the spins of its nuclei are readily polarized, in a process involving contact with rubidium vapor irradiated with a laser beam. In such "hyperpolarized" xenon, the excess of spin-up nuclei can be as much as 20 percent, which gives a far stronger signal than hydrogen's .001 percent spin-up excess. Hyperpolarized xenon has a much longer relaxation time than hydrogen.

[Forward to a Friend](#)

More Industry Headlines

[Senate Agrees to Restore Doctors' Pay](#)

Senator Kennedy's dramatic return results in key vote in the measure to return physician Medicare payment to pre-cut levels.

[Siemens Job Cuts Set at 16,750](#) Company confirms it will cut about 4% of its workforce worldwide; healthcare division to lose 2,800 positions.

[Florida Event a Major Happening](#) Florida International Medical Exposition (FIME) 2008 will help exhibitors gain a competitive edge in the marketplace; Visit DOTmed at Booth #1511.

[MDS Commences Arbitration Against AECL Over Cancelled MAPLE Project](#) Medical isotope maker files \$1.6 billion suit against AECL and the Government of Canada.

[New Legislation Introduced to Override Landmark Supreme Court Decision on Medical Devices](#) Law would allow private state litigation against manufacturers.

[PHILIPS Intera Achieva 3.0T -- Upgrade to a 3.0T MRI!](#) Also see the other great lasers, imaging systems, and more... all on your favorite website for used medical equipment!

[Cardinal Health Expands Manufacturing Capabilities for PET Tracers](#) Three additional Positron Emission Tomography (PET) radiopharmaceutical manufacturing facilities will be opened this year.

[AdvaMed Applauds Commitment to Improve Japanese Patient Access to U.S. Medical Technology](#) Bush administration initiative hailed.

[New Medicare Rule Ensures Access to Care for Beneficiaries in Rural Areas](#) CMS proposed that Medicare beneficiaries who live in rural areas will be able to get their healthcare services from Rural Health Clinics.

[European Commission Working for Cross-Border Healthcare](#) Plan gives greater freedom in choosing treatment options.

Add the ability to associate a single xenon nucleus with a specific molecular target, for example a protein or sugar on the surface of a cancer cell. To do this, the Pines and Wemmer labs have created biosensors equipped with cages that take up and hold onto xenon atoms; the cages, molecules called cryptophanes, are linked to ligands that target specific molecules of interest. Xenon biosensors engineered with several different ligands can be used at the same time; once in place, biosensors carrying hyperpolarized xenon can localize the MRI signals from a range of different molecules on the target.

The final advance underlying the new technique is called Hyper-CEST: hyperpolarized xenon chemical-exchange saturation transfer. While biosensors can bring the xenon to specific molecular targets, in realistic applications relatively few of these are present, only about one percent compared to the total amount of free xenon injected near that target. The signal from the polarized xenon inside the biosensor cages is consequently much fainter than that from the uncaged polarized xenon nearby.

The trick then is to depolarize the xenon nuclei in the immediate vicinity of the cages, which will serve to outline the target in high contrast against the surrounding hyperpolarized xenon pool. This is done through chemical exchange, as xenon atoms are constantly entering and leaving the biosensor cages.

A polarized xenon atom from the pool enters the cryptophane cage, which alters the xenon's resonance frequency, allowing it to be depolarized by rf radiation tuned to a specific frequency. The depolarized xenon atom is then exchanged for a new, incoming polarized atom and reenters the pool. In this way the buildup of nearby depolarized nuclei quickly outlines the target.

Because it produces a much stronger signal, Hyper-CEST acquires images thousands of times faster than would imaging the caged xenon directly. Yet it retains the great advantages of cryptophane biosensors, including their ability to "multiplex," or detect different targets at the same time.

"Slight differences in cage composition, involving only a carbon atom or two, affect the frequency of the signal from the xenon and produce distinct peaks in the NMR spectrum," says team member Tyler Meldrum, of the Materials Sciences Division. "If we design different cages for different xenon frequencies, we can put them all in at once and, by selectively tuning the rf pulses, see peaks at the frequencies corresponding to each kind of cage."

The processes described above - hyperpolarizing the xenon, caging it in biosensors, and building up depolarized xenon in the immediate vicinity of the target through chemical exchange and selective bursts of rf radiation - led to the development of Hyper-CEST MRI. But until now, Hyper-CEST MRI has only been tested at room temperature.

Using biosensor cages as temperature-controlled molecular depolarization gates makes Hyper-CEST MRI possible at a range of higher-than-room temperatures. Because the technique regulates the exchange rate of hyperpolarized-to-depolarized nuclei through the cages, biosensors regulated this way have been nicknamed "transpletors," by analogy to the transistors that act as gates for the flow of electrons from source to drain in electronic systems.

Hyper-CEST at a range of temperatures has many advantages. Most basic is that biomedical MRI must operate at body temperature. Aside from this practical consideration, temperature determines the rates at which different kinds of cryptophane- cage hosts react with their xenon-atom guests. And increasing temperature dramatically increases chemical exchange rates.

The ability to achieve high-contrast images, multiplexed to identify a range of molecular targets, and to do so in a short time, offers

many benefits to patients and physicians.

"Doctors attempting to characterize tumors very often have to take biopsies, and that's painful for the patient, so they usually prefer to take only one biopsy," says Schroeder. "But then they have to run all their tests on this very little tissue. So they would be happy with a method where you have a toolbox of sensors, you throw them all in and wait to let them bind, and then do your tests at the different frequencies and you see what sensors are present, detecting the different proteins. We showed that the exchange rate is so high at increased temperature that you can use a very selective rf pulse."

Enabling fast, sensitive, molecule-specific NMR and MRI in humans and other living subjects is perhaps the most evident advantage of the new technique, but possible applications don't end there. For example, the method offers a better way to study chemical exchange in nanostructures like zeolites, which are important in catalysis, or in versatile carbon nanotubes. Temperature-controlled depolarization is a breakthrough for NMR and MRI that will find uses in a variety of fields.

This research was supported by the Department of Energy's Office of Science, Office of Basic Energy Sciences; by the Deutsche Forschungsgemeinschaft; and by the University of California's Biotechnology Research and Education Program. Berkeley Lab is a U.S. Department of Energy national laboratory located in Berkeley, California. It conducts unclassified scientific research and is managed by the University of California. Visit our website at <http://www.lbl.gov>.

Please [Send us your Comments.](#)

[DOTmed.com](#) | [News Archives](#) | [Advertise](#) | [Unsubscribe](#) | [Contact Us](#) | [Help](#)

Access and use of this site is subject to the terms and conditions of our [LEGAL NOTICE & PRIVACY NOTICE](#)
Property of and Proprietary to DOTmed.com, Inc.
[Copyright ©2001-2008 DOTmed.com, Inc.](#)
ALL RIGHTS RESERVED